EN Product description

Myriad™ is a collagen matrix with an intact extracellular matrix (ECM). Derived from ovine (sheep) forestomach tissue, this advanced matrix is non-reconstituted collagen, thus it retains the innate biological structure and function of the native ECM associated macromolecules including elastin, fibronectin, glycosaminoglycans and laminin.

Myriad™ transforms into a soft conforming sheet, which is naturally incorporated into the surrounding tissue over time. The device is supplied sterile and dry in a variety of sizes and thicknesses, which can be trimmed to meet the individual patient's needs.

Not made with natural rubber latex.

Intended use

Myriad™ is intended for applications in plastic and reconstructive surgery or to cover, protect, and provide a moist wound environment. The device may be fixed, via sutures, staples, or tacks to the surrounding tissue, if desired.

Indications for use

Myriad™ is indicated for implantation to reinforce soft tissue where weakness exists in patients requiring soft tissue repair or reinforcement in plastic and reconstructive surgery, or for management of the following wounds:

- Partial and full thickness wounds
- Pressure ulcers
- Venous ulcers
- Diabetic ulcers
- Chronic vascular ulcers
- Tunneled/undermined wounds
- Surgical wounds (donor sites, grafts, post Moh's surgery, post laser surgery, podiatric, wound dehiscence)
- Trauma wounds (abrasions, lacerations, second-degree burns, and skin tears)
- Draining wounds

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Contraindications

- **Myriad™** is derived from an ovine (sheep) source and should not be used on patients with known sensitivity to ovine (sheep) derived material.
- **Myriad™** is not indicated for use on third degree burns.

Precautions

- **Myriad™** may not have sufficient strength to support stresses encountered in some ventral hernias or large-area, body-wall repairs. **Myriad™** is not intended to augment breast reconstruction.
- Do not apply **Myriad™** in the presence of uncontrolled clinical infection, acute inflammation, excessive exudate, or bleeding.
- **Myriad™** is supplied sterile. Do not use if the pouch seal is broken.
- Discard the device if mishandling has caused possible damage or contamination.
- Single use product. Do not attempt to re-sterilize. Discard all unused portions. Do not reuse.
- Minimize manipulation of **Myriad™** during rehydration and placement.
- Always handle **Myriad™** using aseptic technique.

Potential Complications

The following complications are possible. If any of the complications occur and cannot be resolved, consider removal of the device:

- Infection
- Acute or chronic inflammation
- Allergic reaction
- Seroma formation

Storage

Myriad™ should be stored at room temperature in a clean and dry area.

Directions for use

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

- 1) Rehydrate the device** in sterile saline or sterile Lactated Ringer's solution for a minimum of 5 minutes.
- 2) Trim the device** to fit, if necessary, providing an allowance for overlap. Position the device to achieve maximum contact between the device and surrounding tissue. To facilitate cell migration and tissue ingrowth, an overlap of 3-5 cm with healthy well-vascularized tissue is suggested.
- 3) If desired, suture, staple, or tack into place**, avoiding excess tension. Alternately, the device can be secured using an appropriate secondary dressing.
- 4) Discard any unused portions** according to institutional guidelines for medical waste.

Symbols

ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements		
	5.4.2	Do not reuse
	5.2.6	Do not re-sterilize
	5.2.8	Do not use if package is damaged
	5.1.4	Use by expiration date YYYY/MM/DD
	5.3.4	Keep dry
	5.2.3	Sterilized with ethylene oxide
	5.4.3	Consult instructions for use
	-	Number of Units

AROABIOSURGERY LTD,
2 Kingsford Smith Place, Airport Oaks,
Auckland 2022, New Zealand
Ph: +64 9 869 3035
www.aroabio.com

Endoform® is a registered trademark of **Aroa Biosurgery Ltd.**
Myriad™ is a trademark of **Aroa Biosurgery Ltd.**

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©2019 Aroa Biosurgery Ltd